

## MS-200 MAGNETIC STIRRER

The MS-200 magnetic stirrer is designed to stir liquids placed in vessels on it. It can stir liquids with volume up to 2000 ml.

The stirring speed may be adjusted up to 1250 rpm.

User may set the requested stirring time from 0 to 999 minutes.



### TECHNICAL DATA

Maximal stirring volume:	2000 ml
Stirring speed:	up to 1250 rpm
Top plate size:	135 x 135 mm
Top plate material:	stainless steel
Timer:	0 to 999 minutes
Power requirements:	230V/50Hz
Dimensions:	230 (L) x 180 (W) x 120 (H) mm
Weight:	2.2 kg

**ELMETRON® Sp. j.**

41-814 Zabrze . Witosza 10 POLAND

tel. +48 32 2738106 fax +48 32 2738114

[www.elmetron.pl](http://www.elmetron.pl) e-mail: [info@elmetron.com.pl](mailto:info@elmetron.com.pl)